

# The Dermillume Pro Hand Held IPL!

Red and Blue for Acne Treatments, Red for Rejuvenation Services

## Improve Acne Treatments with affordable phototherapy

- ☞ **Blue Light** kills bacteria deep in the skin where topical treatments can't reach
- ☞ **Red light** promotes healing of lesions by increasing the number of new skin cells in the damaged area

Activates generic ALA  
(aminolevulonic acid) in 45 sec.



Dermillume Pro HH IPL in the Blue-only program

---

## Offer both Acne Treatments and Anti-aging Treatments with One Device

- ☞ Use in **Red + Blue** mode for a *drug-free, effective Acne treatment*
  - ☞ Use in **Blue-only** mode with *generic ALA* (aminolevulonic acid) in dermatological settings
  - ☞ Use in **Red-only** mode for *rejuvenation treatments*
  - ☞ Use in **Red-only** mode to *enhance collagen treatments* – visible results in one week
- 



Dermillume Pro HH IPL in the Red-only program



## The Dermillume Advantage

- ☞ **60 Joules/cm<sup>2</sup> at skin in 90 seconds**
- ☞ **Treats 2.5" x 3" area**
- ☞ **Three programs deliver Blue-only, Red only or Red + Blue modes**
- ☞ **Use to activate Levulan™**
- ☞ **Does not touch skin**
- ☞ **Easy to use – minimal training**
- ☞ **Cordless Design**

Care Electronics, Inc  
303-444-2273 or 888-444-8287

## Optical Specifications for the Pro Hand Held IPL™:

Blue wavelength – 414 nm

Red wavelength – 660 nm

- ☞ **Delivers 9.0 Joules/ cm<sup>2</sup> Blue-only at the skin in 45 seconds for ALA activation**
- ☞ **Delivers 42.0 Joules/ cm<sup>2</sup> Red-only at the skin in 90 seconds for Rejuvenation treatments**
- ☞ **Delivers 60.0 Joules/ cm<sup>2</sup> Red + Blue at the skin in 90 seconds for drug-free acne treatments**

More information at [www.dermillume.com](http://www.dermillume.com)

Levulan is a trademarked product  
of Dusa Pharmaceuticals